

Passivated Tracer Particles with PLL-g-PEG

Developed by: Tim Sanchez

Updated January 2016 by: Stephen J. DeCamp

This protocol passivates 200nm fluorescent sulfate beads and 3um sulfate beads. They are observed to be 100% non-sticky and freely diffusing when mixed in 2mM MgCl₂ buffer (usually makes normal beads sticky), 0.8% PEG, and 1mg/ml microtubules.

PEG was from Laysan Bio: mPEG-SCM, MW 20kDa

PLL was Sigma: number P7890 MW, 15-30kDa

Prepare PLL-grafted-PEG Stock

Want molar ratio of approximately 3.5 Lysine : 1 PEG
Lysine: MW 137Da

- 1) Make 50mM sodium borate buffer, pH 8.5
 - a. 1.906g into 100ml
- 2) Dissolve 1.7mg of p-lysine and 71.4mg PEG into 1.25ml of sodium borate buffer
 - a. Final concentrations: 10mM Lysine and 2.86mM PEG
- 3) Fill up small dialysis membrane to capacity (so that volume does not change much).
 - a. Use 2x 500uL capsules with pore size 5kDa.
 - b. Dialyze against DI H₂O overnight, then replace H₂O, let go another night.
- 4) Recover and aliquot.
 - a. For long term storage, flash freeze and store in -20C.

Coating microspheres with PLL-PEG

The PLL is positively charged, so you must use beads with a negative surface charge.

Carboxylated beads supposedly work, according to this protocol:

(Title: **Ligand-Specific Targeting of Microspheres to Phagocytes by Surface Modification with Poly(L-Lysine)-Grafted Poly(Ethylene Glycol) Conjugate**).

Sulfate beads are supposed to be very negatively charged, and thus work very well.

Ratio of PLL-PEG:beads should be dependent on the surface area of the beads.

3µm beads:

- 1) Mix 10µL PLL-PEG(~1.37mg/ml of PLL) + 100µL HEPES(100mM) + 12µL beads (10%) + 878µL H₂O
- 2) Incubate an hour on a rotating platform.
 - a. **Ultra-sonicating at 20% (low setting) inside the 1.5ml eppendorf tube was very helpful for the 3µm beads and ESSENTIAL for the 200nm beads. This prevents aggregation into clusters of beads.**
 - b. Sonicate once every 20 minutes for the 1 hour.
 - i. Each time, sonicate 3X for 3 second bursts (too much can melt the beads).
 - c. Put back on rotating platform.
- 3) Transfer to 500µL eppendorf tubes. Spin at 4k (5140g) for 5 minutes. Will get some beads glazing the walls, but should get a reasonable pellet, too.
- 4) Aspirate supernatant and consolidate pellets to another 500µL tube.
- 5) Exchange buffer (whatever buffer you like). Vortex and pellet should resuspend easily.
- 6) Spin again as in step 3.
- 7) Aspirate and resuspend in whatever volume you want. (approximately 50µl)
 - a. Transfer to new tube for storage in 4C.

200nm beads:

- 1) Mix 15µL PLL-PEG + 20µL HEPES (100mM) + 2.4µL bead stock + 162.6 H₂O
- 2) Incubate an hour on a rotating platform.
 - a. **Ultra-sonicating at 20% (low setting) inside the 1.5ml eppendorf tube was very helpful for the 3µm beads and ESSENTIAL for the 200nm beads. This prevents aggregation into clusters of beads.**
 - b. Sonicate once every 20 minutes for the 1 hour.
 - i. Each time, sonicate 3X for 3 second bursts (too much can melt the beads).
 - c. Put back on rotating platform.
- 3) Transfer to 500µL eppendorf tubes. Spin at 14k (18,000g) for 5 minutes. Will get some beads glazing the walls, but should get a reasonable pellet, too.
- 4) Aspirate supernatant and consolidate pellets to another 500µL tube.
- 5) Exchange buffer (whatever buffer you like). Vortex and pellet should resuspend easily.
- 6) Spin again as in step 3.
- 7) Aspirate and resuspend in whatever volume you want. (approximately 40µl)
 - a. Transfer to new tube for storage in 4C.